•	Chapter 7:
•	程式设定   Program (PT) Test

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
# Chapter 7 程式设定 | Program (PT) Test

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<ul style="list-style-type: none"> <li>④</li> </ul>	湿度区域 Humidity	湿度待机区 Humidity standby area
<ul style="list-style-type: none"> <li>Waiting area: if it is temperature, it refers to the temperature area, and if it is humidity, it refers to the humidity area. <a href="https://www.wewontech.com/voc-emission-chamber/">https://www.wewontech.com/voc-emission-chamber/</a></li> </ul>		The Relationship Between Standby Action and Standby Time:
1.1 The actual value within the standby time reaches the standby range, the release of standby action		
1.2 The actual value within the standby time does not reach the standby range, the release of standby action		

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			program, set the block where the partial cycle starts running.
•	③	开始段号   Start (Segment)	已设置程序中，设置部分段循环运行开始的程序段   Set the number of cycles of the program. When it is 0, it will loop indefinitely.
•	④	结束段号   End (Segment)	已设置程序中，设置部分段 循环运行结束的程序段，小于 0 时不循环   In the set program, set the section of the loop that ends the loop. If it is less than 0, it will not loop.
•	⑤	循环次数   Number(Recycle Times)	已设置程序中，设置部分段 循环运行的循环次数，小于 0 时不循环   In the set program, set the number of cycles in which the partial segment is cyclically run. When it is less than 0, it does not cycle.
•	⑥	参数确定   Save	把当前参数输入控制器   Input current parameters into the controller
•	⑦	参数上传   Upload	把当前参数上传显示器   Upload the current parameters to the display
•	⑧	连接到   Link	当前程序运行结束后要连 续运行程序的编号   The number of the program to be continuously run after the current program finishes running.

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设定   SETTING		讯控动作   Signal Control Action
SEG TIME ≥ TIME ON + TIME OFF	= 0.00 Signal Control # SET 4	<p>ON</p> <p>Signal 1</p> <p>OFF</p> <p>Time Off</p> <p>程序段 SEG1 SEG2 SEG3</p>
	TIME ON = 1.00 Signal Control # SET 5	<p>ON</p> <p>Signal 1</p> <p>OFF</p> <p>Time On</p> <p>Time Off</p> <p>程序段 SEG1 SEG2 SEG3</p>
SEG TIME < TIME ON + TIME OFF	TIME ON = 3.00 Signal Control # SET 6	<p>ON</p> <p>Signal 1</p> <p>OFF</p> <p>Time On</p> <p>Time Off</p> <p>程序段 SEG1 SEG2 SEG3</p>
	TIME ON = 3.00 Signal Control # SET 7	<p>ON</p> <p>Signal 1</p> <p>OFF</p> <p>Time On</p> <p>Time Off</p> <p>程序段 SEG1 SEG2 SEG3</p>